NEW PORTABLE ISOLATION ROOM FOR HOSPITALS BY IBIOM

The PIRH can be universally attached to most transport stretchers on the market.

Used to help reduce the spread of the virus when a patient is transported to dedicated units.



- Once at the destination, the PIRH is removed and cleaned so the patient can be transferred.
- The PIRH can be easily removed by one or two caregivers due to its rigid structure and its light weight of 17 lb (7.8 kg).
- The tubing can be passed under the vinyl to provide IV fluids or ventilation to the patient.
- PIRH can be stacked for easy storage.

The PIRH has an inside height of 31.75 in (80.6 cm) allowing a clearance of approximately 26.5 in (67.3 cm) inside for the patient.

The width at mattress level is 26 in (66 cm) suitable for stretchers with a 26 in (66 cm) mattress. A 30 in (76 cm) version is also available.

- No ventilation or negative pressure is included.
- Never leave a patient unattended under the PIRH.
- A ventilator may be required depending on the patient's condition.

