

HEMA+®

Treatment chair & phlebotomy chair

HEMA+ Must configuration

Included features:

- » 2 functions hand control
- » Standard column version: 43 to 68.5 cm / 17" to 27"
- » 5.1 cm / 2" Seamless cushions
- » Electric Trendelenburg -5°
- » 12.7 cm / 5" Twin casters with central lock

HEMA+ Comfort configuration

Included features:

- » HEMA+ Must configuration
- » Rechargeable battery pack*
- » Manual Trendelenburg handle (-5°)*

Armrest options

Must select one option, always in pairs:

- » Padded, height adjustable and swivel armrests (15.2 x 38.1 cm / 6" x 15")
- » Retractable armrests (10.1 x 25.4 cm / 4" x 10")
- » Non adjustable armrests with retractable shelves (15.2 x 30.5 cm / 6" x 12")

Additional option available

- » High column version: 51 to 91.5 cm / 20" to 36"

Accessory

- » 325M: Silicone armrest

Specifications

Umano Medical reserves the right to change specifications without notice.

Model

- HEMA+

Adjustments

- Backrest: 0° to 80°
- Legs: 0° to -100°
- Armrests:
 - Non adjustable with retractable shelves
 - Height adjustable and swivel
 - Retractable
- Trendelenburg: -5°

Maximum weight capacity

- 159 kg / 350 lb using medical-grade electrical components
- 120V - 3.9 AMP - 60 Hz

Dimensions

- Standard width: 50.8 x 162.6 cm / 20" x 64"
- Backrest: 80 x 50.8 cm / 31.5" x 20"
- Seat: 45.7 x 50.8 cm / 18" x 20"
- Legs: 30.5 x 50.8 cm / 12" x 20"

Cushions

- 5.1 cm / 2" Thick
- Seamless for easy infection control
- Vinyl upholstery resistant to common cleaning products¹

Controls

- 2 functions hand control

Casters

- 12.7 cm / 5" Twin casters with central lock

* Optional

1. Internal data on file



HEMA+

Treatment chair & phlebotomy chair



Umano Medical Inc.
230, boulevard Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll Free: 1 (844) 409-4030
information@umanomedical.com
umanomedical.com

© Umano Medical Inc. 2026. All rights reserved - 905-9000-ENG R2 - Printed in Canada
The Umano Medical and HEMA+ marks and/or logos are trademarks owned by Umano Medical World Inc. and used with permission by Umano Medical Inc.

HEMA+®

Treatment chair & phlebotomy chair

In use: phlebotomy, nuclear medecine, ophtalmology, emergency and psychiatry. Designed for ergonomic positioning during a variety of short-duration treatments or procedures, including drawing blood and general patient exams. Easily accessible to patients of all sizes. Prevents patient falls.

Cleaning: The vinyl upholstery is tested for common cleaning products used in hospitals and clinics, which is imperative when infection control is a priority¹.



- Colors:
- Denim Blue
 - Raven
 - Other on request


Entirely designed and manufactured in Canada, locally manufactures offers you easy access to spare parts and customer service in English and French.



Silicone armrest.



2 functions hand control.



Quick-release manual handle.
To bring patient into Trendelenburg position.